



# SWISS BIMANUAL I/A



# SWISS BIMANUAL I/A



## The best bimanual choice

The SWISS BIMANUAL I/A handpieces are characterised by their well-established aeroplane shape. The textured handpiece grip surface offers comfortable handling and user control. They are compatible with all common luer connectors and available in the preferred designs and sizes 23G and 21G.

**SWISS BIMANUAL I/A systems (single-use / sterile)** consisting of: Bimanual I/A instrument (box of 10)

Item Code	Size	Incision Size*	Aspiration Cannula	Aspiration Port	Aspiration Head	Irrigation Cannula	Irrigation Port
SBS310A	21G / 0.8mm	1.1mm	aeroplane shape	0.35mm	textured	aeroplane shape	2 ports, 0.45mm
SBS320A	23G / 0.65mm	0.8mm	aeroplane shape	0.3mm	textured	aeroplane shape	2 ports, 0.35mm

aeroplane  
shape




**SWISS BIMANUAL I/A systems (single-use / sterile)** consisting of: Bimanual I/A instrument (box of 10)


Item Code	Size	Incision Size*	Aspiration Cannula	Aspiration Port	Aspiration Head	Irrigation Cannula	Irrigation Port
SBS105	21G / 0.8mm	1.1mm	curved aeroplane shape	0.35mm	textured	open	1 port, 0.5mm
SBS110A	21G / 0.8mm	1.1mm	curved aeroplane shape	0.35mm	textured	curved aeroplane shape	2 ports, 0.45mm
SBS123A	23G / 0.65mm	0.8mm	aeroplane shape	0.3mm	textured	curved aeroplane shape	2 ports, 0.35mm

\* Incision sizes are indications only based on average experiences. The individually required incision sizes may vary depending on surgical techniques

SBS105



SBS110A  
SBS123A



instruction for use



medical

Medicel AG, Dornierstrasse 11, 9423 Altenrhein, +41 71 727 10 50, medical.com



digital brochure

made in switzerland