


ophtec

IOLMATIC®

Injector System for
Precizon™ IOLs



Ready to use with
just 2 simple clicks!



Injector System for Precizon™ IOLs

IOLMATIC®

Our new ready to use packaging for Precizon™ IOLs, is fully compatible with the IOLMATIC® injector system. **With just 2 simple clicks** it is connected to the lens holder, loaded with the IOL and ready to use. It's easy, saves time and ensures a safe IOL delivery with reliable results.



View full
instruction
video

Get started

Your sales representative is happy to provide a short training for using the IOLMatic® Injector System. For more information and the instruction video, please go to ophtec.com/iolmatic



Specifications

- For single use only
- 2.4 mm incision
- Push injection system
- Compatible with Precizon™ IOLs
- Soft tip for 1 piece acrylic IOL
- Product code:
IOLMATIC 2.0-A (1 unit)
IOLMATIC injector system 2.0 S (1x1)